

PRODUCT UPDATE

Motion Control
division of
Fillauer

115 N. Wright Brothers Dr. • Salt Lake City, UT 84116 • USA
888.696.2767 • 801.326.3434 • 801.978.0848 (f) • info@UtahArm.com • www.UtahArm.com

Spring 2016

“ETD2” Preview



- Shorter than original ETD by 3 cm
- Secure gripping of flat and cylindrical objects
- FLAG option and Bluetooth®
- Auto-Cal - triggered by iPhone®
- Outside hook friction pads
- Replaceable finger lining
- Multi-Flex or Flexion Wrist options
- Safety release

Available Summer, 2016
(limited release)

MC ProPlus Hand

- FLAG option and Bluetooth®
- Heavy-duty, brushless, DC motor
- Regal silicone glove option
- Multi-Flex or Flexion Wrist option
- Safety release
- Auto-Cal - triggered by iPhone®
- Short hand option



FLAG

FORCE LIMITING AUTO GRASP

FLAG is a simple, but very useful feature, now available in ProPlus versions of MC Electric Terminal Device (ETD) and MC Hand.

- **Force Limiting** - Grasp a delicate object and the terminal device will stop gripping at ~ 2 lb of gripping force; each 'pulse' raises gripping force by ~ 2 lb
- **Auto Grasp** - Auto Grasp detects socket or electrode 'slip', and automatically increases grip force

Suggested L-Codes: L6881 (Auto Grasp); L6882 (Microprocessor Control)

Wireless Bluetooth®



For iPhone® 4S and above, iPad® 3rd gen. and above, iPod touch® 5th gen. and above

- Connect to Apple® hand held devices
- The iOS App (MCUI) is **FREE** for your Apple® device at the App Store
- Auto-Cal (wearer-triggered auto adjustment)

COAPT Pattern Recognition Now Available from Motion Control

COAPT has now named Motion Control an official distributor of COAPT Products.

Pattern Recognition represents a new era in control systems for myoelectric devices. With the addition of Pattern Recognition, operating the prosthesis becomes intuitive and natural. Switching between degrees of freedom is virtually seamless.



Motion Control Products compatible with COAPT

- Utah Arm U3+
- Utah Hybrid Arm
- MC Standard Wrist
- MC ProWrist
- MC Terminal Devices

Our Online Catalog is LIVE!

We've organized all our products in a convenient, online catalog.

Our Products are now accessible from your:

- Desktop or laptop computer
- Hand-held tablet
- Smart phone

In fact, anywhere the internet is available you can see, and show your patients, a photograph and description of our products.

Features:

- Description of each part
- Photograph of each part
- Mouse-over zoom for a closeup view of connectors and adapters (see inset)
- Place your order online (not available in all countries); or simply submit the order to us for advice
- U.S. List Price with Log In (for Authorized Registered Users)



www.UtahArm.com/Catalog

Search by:

- Category
- Keywords
- Part Number



115 N. Wright Brothers Dr. • Salt Lake City, UT 84116 • USA
888.696.2767 • 801.326.3434 • 801.978.0848 (f) • info@UtahArm.com • www.UtahArm.com