



# UTAH ARM 3+

**N**ew and exciting features have been Added to the Utah Arm, evolving the



Our goal is easier usage for patient and prosthetist, with the new Dual Lock System, and Wireless Connection to the computer.

## **NEW!**

- Dual Lock System: friction lock plus the Hi-Load Lock Pin
- BlueTooth Connection
- Silent FreeSwing
- Functional MC TD's with Flexion, Multi-Flex, In-Hand Wrist Rotator
- iLimb Hand Compatible

## **USUAL HIGH SERVICE:**

- Versatile Input Choices
- Rapid Service & Dependability
- Borrow a Demo or Loaner
- Price - Still affordable

*Motion Control*  
KEEPING LIFE IN MOTION

Call 1.888.MYO.ARMS or visit [www.UtahArm.com](http://www.UtahArm.com)



# Adds a New Level of Performance:

Impressive new features of the Utah Arm 3+ (U3+) will add ease of use — and ease of fitting for the prosthetist.

Look to the U3+ for *function, ease of use, affordability, and service!*

### CHECK IT OUT:

### WHAT IT OFFERS:

• Dual Lock System...	Elbow stops momentarily – brake engages Hand function is immediately available – No Wait! High load on elbow (~10lb) engages the lock pin
• Wireless Connection...	BlueTooth with Hard-Wire Backup No modification to U3+ Arm – Plug & Play!
• Silent FreeSwing...	Totally silent – no power consumed! Automatic at full extension
• iLimb Hand Compatible...	Also – MC TDs with Multi-Flex Wrist, Flexion and high-speed MC Wrist Rotator Virtually all TDs on the market compatible with U3+

### UNMATCHED MOTION CONTROL SERVICE:

• Versatile Input Choices...	New Triad Preamps Use snap electrodes when needed Harness-mounted sensors provide convenience
• Rapid Service & Dependability...	Call Us – WE'RE LIVE! Rapid turnaround on orders & repairs
• Borrow a Demo / Backup...	Educate your patient and clinic Evaluate all options Rental Plan – Unique Trial Fit Opportunity
• Price – Still affordable...	Only a modest premium for New U3+ features

